

# Project G5: Melanoma Detection using Knowledge Distillation for Mobile Phones

- **Validate teacher complexity:** ResNet-34 vs ResNet-50 vs EfficientNet-B2 on validation accuracy first.
- **Simplify hyperparameter search:** Proposed grid ( $T=[1,2]$ ,  $\alpha=[0.1,0.5,0.7,0.9]$ ,  $LR=[0.00001-0.01]$ ,  $batch=[16-256]$ ,  $dropout=[0.1-0.5]$ ) = hundreds of configs—**use Bayesian optimization we discussed in class** or prioritize  $T$  and  $\alpha$  only.
- KD is appropriate for mobile deployment, but **establish teacher performance ceiling** (train ResNet50, measure AUC/ECE) before distillation—if teacher  $ECE > 0.1$ , fix calibration first.

## Actions before Milestone-2

- **Dataset splits:** Download HAM10000, create 70/15/15 train/val/test (stratified by melanoma label), **not 80/20 train/test with K-fold CV**—CV is too expensive for KD experiments.
- **Simple baselines:** Train logistic regression on metadata (age, sex, lesion location) + handcrafted features (color variance, asymmetry), expect  $\sim 70\%$  AUC baseline.
- **Train teacher first:** ResNet50 (ImageNet pretrained, frozen backbone, fine-tune head) with weighted BCE (weight=8 for melanoma), train  $> 20$  epochs, log AUC, ECE, and precision@95% sensitivity on val.
- **Temperature scaling:** If teacher  $ECE > 0.05$ , apply post-hoc temperature scaling on validation set before KD—properly calibrated teacher is essential for distillation.
- **Recommended tracking setup:** W&B with fixed seeds, log loss curves (KL + BCE components separately), model size (MB at fp32/fp16/int8), inference latency (ms on CPU), memory footprint (MB RAM).

## Ablations

- **Temperature  $T$ :** Compare  $T=1$  (hard targets) vs  $T=2$  (softer targets)—hypothesis is  $T=2$  improves student calibration (lower ECE by 0.02–0.03) following Islam et al.
- **Loss weighting  $\alpha$ :** Test  $\alpha=0.5$  (balanced) vs  $\alpha=0.9$  (teacher-heavy)—expect  $\alpha=0.5$  yields best AUC-ECE tradeoff, avoiding over-reliance on uncalibrated teacher.
- **Quantization:** Post-training quantization (fp32→int8) should drop AUC by  $< 2\%$  and ECE by  $< 0.03$ —if degradation higher, try quantization-aware training (QAT).

## Risks & mitigations

- **Risk:** HAM10000 has 11% melanoma imbalance, student may underfit minority class—**Mitigation:** Use weighted BCE (weight=8 for mel), focal loss ( $\gamma=2$ ), or oversample melanoma 2:1 during training.
- **Risk:** Extensive hyperparameter grid ( $T$ ,  $\alpha$ , LR, batch, dropout, quantization) is computationally prohibitive—**Mitigation:** Prioritize  $T$  and  $\alpha$  (4 configs), fix other hyperparams to literature defaults (LR=0.001, batch=64, dropout=0.2).
- **Risk:** Quantized student loses calibration (ECE spikes)—**Mitigation:** Measure ECE on quantized model separately, retrain with QAT if post-training quantization  $ECE > 0.1$ .

## Open questions

- How will you handle domain shift from HAM10000 (dermoscopy, clinical lighting) to user-taken smartphone photos (variable lighting, distance, angle)—will you test on any real-world images?
- Will you ablate global max pooling vs global average pooling in teacher ResNet50 before distillation, following Islam et al.’s finding that max pooling outperforms average?
- What specific quantization method (PyTorch native, TensorFlow Lite, CoreML)—and will you measure calibration (ECE) separately on quantized student vs fp32 student?